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Capacitive Resistive Monopolar Radiofrequency (CRMRF) therapy at 448 kHz: The effects on deep blood flow and elasticity of tissues

Binoy Kumaran, Anthony Herbrand & Tim Watson, School of Health and Social Work, University of Hertfordshire, United Kingdom

www.herts.ac.uk b.r.kumaran@herts.ac.uk

Introduction

Radiofrequency-based electrophysical agents (EPAs) have been employed in therapy practice for several decades. They are used to reduce pain and inflammation and enhance tissue healing. Besides **Shortwave Therapy (SWT)** that operates at **27.12 MHz**, which is the commonest RF-based EPA, devices operating at significantly lower frequencies are also available claiming to be clinically effective although their evidence is limited.

CRMRF is one such EPA that operates at a base frequency of **448 kHz** (0.45 MHz) and delivers RF in continuous mode.

Purpose

This study aimed to investigate the changes in **blood flow** and **elasticity** in tissues at least at two centimetres depth from the skin in response to cutaneous application of CRMRF (administered using **Indiba Activ 902**; Indiba S. A., Barcelona, Spain) in healthy adult volunteers, in order to examine its physiological effects.

Methods

Participants: Seventeen healthy volunteers (10 females, 7 males) randomly selected from among the members of the University of Hertfordshire took part in this study.

Mean ($\pm$ SD) subject characteristics

Age (years)	Height (m)	Weight (kg)	Body fat (%)	Visceral fat	Body Mass Index
45.7 (12.7)	1.7 (0.1)	71.5 (10.0)	30.3 (7.6)	7.2 (2.5)	24.7 (2.7)


SD – Standard Deviation; kg – Kilogram; m – Metre

Design: Single-blind randomized crossover study comprising four groups, with all 17 participants attending all four groups in a random order – high, low and placebo dose groups and a control group.

Intervention: The high, low and placebo dose groups received CRMRF therapy for 15 minutes given to the medial aspect of the lower end of thigh. The control group did not receive any intervention and simply rested for 15 minutes.


Mean ($\pm$ SD) power delivered in Watts

	High Dose	Low Dose	Placebo Dose
CAP Mode Power	24.0 (3.6)	7.2 (2.3)	2.1 (1.4)
RES Mode Power	51.6 (6.2)	24.6 (5.0)	3.2 (1.4)
Overall Power	42.4 (4.6)	18.8 (3.8)	2.8 (1.2)


Measurements: Baseline and post treatment blood flow velocity, blood volume and tissue elasticity were recorded using Doppler Ultrasound and Ultrasound Elastography (Esaote MyLab70 XVG, Italy). The core temperature, blood pressure and pulse rate were concurrently monitored.

Colour Doppler image showing blood velocity

Power Doppler image showing blood volume



Elastography image showing information on tissue elasticity



Data analysis: Computational analysis of ultrasound images was performed using MATLAB algorithms (MathWorks, Massachusetts), which converted the colour image data into numerical data. Group data were compared using a **two-way repeated measures ANOVA** model using SPSS Statistics (v20). Statistical significance was set at $p \leq 0.05$ (0.8 P, 95% CI).

Results: A significant overall difference in both **volume** ($p=0.047$, 0.82 P) and **intensity** ($p=0.009$, 0.98 P) of blood flow was demonstrated between the four conditions.

No significant differences were noted between the high and low dose groups although the high dose response was greater.

No significant changes were noted in the blood flow velocity measurements or tissue elasticity measurements.

The core temperature, blood pressure or pulse rate did not change significantly at any time.


Conclusions: Both the high and low doses of CRMRF at 448 kHz are capable of significantly increasing the volume and intensity of blood flow at depth, but not the velocity of flow or elasticity of tissues.

References

- Allan PL. Clinical Doppler ultrasound. London: Churchill Livingstone; 2000.
- Al-Mandeel MM, Watson T. The thermal and nonthermal effects of high and low doses of pulsed short wave therapy (PSWT). Physiother Res Int. 2010;15(4):199-211.
- Kumaran B, Watson T. Radiofrequency-based treatment in therapy-related clinical practice – a narrative review. Part I: acute conditions. Phys Ther Rev. 2015;20(4):241-54.

Acknowledgements

The Health and Human Sciences Ethics Committee with Delegated Authority (HHSECD) of the University of Hertfordshire, approved this study (hhsecd@herts.ac.uk). This study was funded by Indiba S. A., Barcelona, Spain. The industry funders had no role in the study design, data collection, data analysis or the preparation of this abstract. Physiotherapy UK 2015 Conference, 16–17 October 2015, Liverpool, United Kingdom.


Clinical research

Adipose- Derived Stem Cell Proliferation (in vitro)

Treatment with 448kHz electric stimulus increased cartilage-specific collagen type II by 51% and glycosaminoglycans by 20% in respect to sham groups.

Hernández-Bule ML, María Trillo Á, Martínez-García MÁ, Abilaloud C, Úbeda A (2017) Chondrogenic Differentiation of Adipose-Derived Stem Cells by Radiofrequency Electric Stimulation. J Stem Cell Res Ther 7: 407


Clinical research

Improving Pain and Function in Osteoarthritis RCT

4 weeks INDIBA treatment with exercise and advice (control and placebo groups) with 3 months follow up

INDIBA treated group had a 66% reduction in pain from baseline, at 3 months follow up pain was still reduced by 45%. Control group was 8%

Functional quality of life was measured by WOMAC; 45% decrease in score at 4 weeks and 38% at 3 months.

[Read more](#)

Kumaran B, Watson T. Treatment using 448kHz capacitive resistive monopolar radiofrequency improves pain and function in patients with osteoarthritis of the knee joint: a randomised controlled trial. *Physiotherapy* (2017).

Clinical research

Blood flow response to INDIBA applications

Skin flow response from 17 participants showed a significant change in the group treated with INDIBA in respect to placebo and control groups.

Kumaran B, Watson T. Capacitive resistive monopolar radiofrequency (CRMRF) therapy at 448kHz: localised application significantly enhances and sustains skin physiological responses. World Conference for Physical Therapy congress 2015, Singapore; 05/2015.